

Instructions for Installation

Operating System: Windows 7

Imaging Program: Xray Vision version 4 (and rebrands)

Note: *Windows 7 is no longer a supported Operating System from Microsoft. Support for Windows 7 was discontinued in 2020 and receives no further security updates to the operating system. This in turn, makes Windows 7 a security risk, and is no longer considered to be HIPPA compliant.*

Files needed: (files can be found at www.tuxedoimaging.com under Support & Downloads)

- Tuxedo A Series Driver
- XV 4

DRIVERS


- [Tuxedo A Series Driver](#)
- [Windows 10 Patch Installer](#)
- [Windows Secure Boot Patch](#)

PLUG-INS

- [XVCapture / DCV 3.5](#)
- [XVCapture / DCV 4](#)
- [XV.4](#)


Step 1: Ensure that the Tuxedo A Series sensor is **NOT** plugged into the PC.

Step 2: Install the Tuxedo A Series Driver (*filename seen below*)



This process may prompt for a restart of the PC, if that happens, restart the PC before moving forward.


Step 3: Plug in the Tuxedo A Series sensor into the computer's USB port directly. You will likely be prompted with a pop-up of Windows searching for the driver software.




- You will need to wait for both green checkmarks to appear, before you can close this prompt. This process may take several minutes.

Step 4: Run Xray Vision (or rebrand) as administrator, let the program open, and then close the program with the 'X' in the upper right-hand corner.

Step 5: Install the Tuxedo Upgrader, by running it as Administrator. (*filename seen below*)



- Click "Next" in the installer window until, and "Yes to All" until "Finish" is no longer greyed out.
- Verify there were no errors.



- Click "Finish"

Step 6: Run Xray Vision (or rebrand) and go into the Toolbar at the top, select **Tools > Hardware Options > TuxedoASeries**, and then use the following recommended setup for each of the tabs in the following screenshots.

General | GBC | Default Calibration | Sensor Configuration | Custom Path for Calibration Files

Show switch sensor button
 Close hardware between each capture
 Enable undereposure detection
 Make comers black Added comer: 0

Enable logging
 Save diagnostic images
 Close hardware after layout capture is completed

What type of image should be returned?
 Low Contrast High Contrast

What type of image processing should be applied?
 Normalize Adaptive Normalize
 Laplace Edge Enhancement Mask Size: 3x3 Factor: 3

Apteryx General Enhancement (Almg.DLL) filters
 Enhance Local Contrast (CLAHE) Scale = 5; Clip = 3.80 Edit
 Adaptive Normalize Low = 0.002; High = 0.002 Edit
 Median Blur Mask Size: 3x3
 Sharpen Mask Mask Size: 19x19 Factor: 15
 Gauss Blur Mask Size: 3x3

General | GBC | Default Calibration | Sensor Configuration | Custom Path for Calibration Files

GBC adjustments are used to compensate for radiation variations in your digital x-ray equipment

Posterior Anterior
 Brightness: 0 Brightness: 0
 Contrast: 0 Contrast: 0
 Gamma Correction: 0.40 Gamma Correction: 0.40

General | GBC | Default Calibration | Sensor Configuration | Custom Path for Calibration Files

Auto rearm Hardware Timeout: 300
 Enable image binning

XRay Detection
 Detection Mode: XVIS detection XVIS detection threshold: 250 mV

Acquisition
 Integration Mode: Max between TWI and X det Integration Time (ms): 3000
 Gain: 1.5

Chain offset
 Chain offset: Computed by COC Manual chain offset (mV): 0

Step 7: Press **OK** until you are back to the main screen.



Step 8: To obtain the calibration files for your Tuxedo sensor, please contact Tuxedo Imaging support.